



Training on LJ culture method

Module 10: Identification And Confirmation Methods

DATE:

VENUE: SRL, Uganda


FACILITATOR:

Content outline

- Definitive *M. tuberculosis* identification
- Cultural tests:
absence of growth on LJ containing PNB
- Biochemical tests:
niacin
nitrate reduction
Catalase

Isolation and identification of mycobacteria

Sample testing flow chat



Isolation


Growth of AFB detected

Identification

- Conventional biochemical tests (4–8 weeks)
- GeneXpert (2 hours)
- Line-probe assay (1-2days)
- **Antigen detection strip (15 minutes)**



M. tuberculosis identification



Phenotypic approach: biochemical tests

- Advantages:
- Inexpensive
- Will identify members of MTB complex
- Gives reliable results for some MOTT.
- Disadvantages:
- Labour-intensive
- Long turn-around times
- Technical expertise required
- Limited to solid cultures

Genotypic approach

- **Advantages:**

1. Fast
2. Able to distinguish between members of MTB complex
3. Identifies NTM
4. Fewer biosafety concerns

- **Disadvantages:**

1. Expensive
2. Requires dedicated equipment
3. Requires technical expertise
4. Traditional methods still required for culture and DST

Phenotypic identification

- Has to consider all the characteristics, including:
- Morphology (colonies/bacilli)
- Culture based
- Biochemical tests

Culture based Technique: Growth on medium containing p-nitrobenzoic acid (PNB)

This is based on Selective inhibition for growth of MTB allowing only NTMs to grow

Procedure

1 LJ slant with PNB at 37 °C

Examine at 28 days

1 LJ slant without PNB at 37 °C

Interpretation


- Abundant growth on both slopes: mycobacterial strain other than tubercle bacilli.
- Abundant growth on control tube and little or no growth on PNB medium: MTB complex strain.
- No growth on either slope: non-interpretable test, to be repeated.




Microscopy Identification

This involves staining with ZN and Examining the morphology of the Bacilli

Cord-like Morphology for MTB



Striped” Morphology for NTM



Immuno- chromatographic assay: Antigen Test


- Immunochromatographic assay: Simple lateral flow speciation test, Monoclonal antibody technology
- Detects MPB64 antigen- predominant protein Ag, secreted specifically by members of the MTBC
- First monoclonal antibody labelled by colloidal gold particles reacts with MPB64 antigen in sample to form antigen-antibody complex
- Complex is then captured by a second monoclonal antibody fixed in the middle of the test zone



Antigen detection using Immunochromatographic Assay

Procedure:


- A 100 µl suspension of growth from liquid or solid media is added to the specimen well.
- The test is read after 15 minutes
- All procedures above should be done in a certified BSC




Antigen detection using Capilia TB

Interpretation:


MTB complex – two purple-reddish bands at test and control sites



NTM or other bacterial species – one purple-reddish band only at control site



Note: The test cannot be interpreted if there is no band at the control site



Advantages of Immunochromatographic Assay

- Performed directly from the AFB positive liquid and solid culture
- Detects MTBC in pure, mixed and contaminated liquid cultures
- Does not require additional equipment or reagents
- Results immediately available

Limitations of Immunochromatographic Assay

- Does not differentiate between members of the MTBC
- Requires culture (no direct inoculation from clinical specimen)
- Some substrains of *M. bovis* BCG are interpreted as negative (organism lacks MPB64)
- Strains of microbes, such as *S. aureus*, that produce protein A may produce a false positive result
- False Negative if the MPB64 concentration is below the detectable limit or if mutations have occurred in the MPB64 gene of MTBC

Exercise (3 minutes)

True” or “false”

1. *M. tuberculosis* can be identified only by the morphology of the colonies.
2. Expired reagents affect the results of identification tests.
3. The Immunochromatographic test allows fast identification of TB complex from positive cultures.

Biochemical tests

These methods are currently less applied, they include

- Niacin production
- Nitrate reduction
- Catalase negative at 68 °C

Niacin test-Principal



This is a Metabolite accumulating in medium during M.tb growth:

- Differentiates M. tuberculosis from other species.
- Rarely positive in mycobacteria other than M.tb.
- The test must be carried out on pure cultures otherwise it will yield false results
- Pure cultures of acid-fast bacilli grown on solid medium, age 21-28 days, with more than 50 colonies

Niacin test - procedure

Commercially available test-
strips used.

Interpretation



Nitrate reductase test

- *M. tuberculosis* reduces nitrates to nitrites.
- Test takes about 4 hours to report results.
- Cultures tested:
 - 4 weeks old
 - abundant growth.

Nitrate reductase test - results and interpretation

- *Negative*: no colour.
- *Positive*: red colour, varying from pink to very deep crimson:

faint pink = +/- (inconclusive)

clear pink = 1+


deep pink = 2+

red = 3+

deep red = 4+

purplish red = 5+

POS




Catalase test: principle

- Catalase converts H_2O_2 to water and O_2 bubbles .
- All mycobacteria produce catalase except for rare isoniazid-resistant mutants of *M. tuberculosis* and *M. bovis*.
- catalase enzyme present in most of the mycobacterial species is heat-stable at pH 7.0 except in tubercle bacilli

Catalase test: Procedure and interpretation

- 1 tube: incubation at 68 °C for 20 minutes- No bubbles hence no MTBC
- 1 tube: incubation at room temperature for 20 minutes- bubbles hence Positive for MTBC



COLOR PLATE 8. Heat stable catalase test.

Assessment

1. State the principal of the Immunochromatographic assay used in MTB identification.
2. List three advantages and two limitations of using Immuno-chromatographic assay test method

Summary

Characteristics of *M. tuberculosis* AFB

- Slow growth rate
- Growth temperature 35-37 °C only
- Typical morphology
- No pigmentation

Additional

No growth on LJ medium containing p-nitrobenzoic acid

Niacin test: positive

Catalase test: negative at 68 °C

Nitrate reductase test: positive



REFERENCES

- www.who.int/tb/laboratory/mycobacteriology-laboratory-manual.pdf
- Grandjean et al. 2008
- Global Tuberculosis Report, WHO 2019
- www.who.int/tb/publications/2012/tb_biosafety/en/
- medicine.kln.ac.lk/depts/publichealth/Fixed_Learning/Campaigns/TB%20Campaign/Manuals/Laborotory/Introduction.pdf
- www.ghdonline.org/uploads/Isolate_storage_packaging_and_transportation
- jcm.asm.org/content/36/2/402
- www.ncbi.nlm.nih.gov/pmc/articles/PMC3838071/

Acknowledgments

