

# **Training on *Mycobacterium tuberculosis* drug susceptibility testing (first and second line LJ DST)**

## **Module 12: Trouble shooting of LJ DST results**

**Venue:**

**Presenter:**

**Date:**

# Introduction and Objective

## Introduction

This module details the steps taken when troubleshooting discordant and invalid/noninterpretable LJ DST results.

## Objective

By the end of the module participants should be able to troubleshoot discordant and non interpretable LJ DST results

# Module Outline

- Use of correct equipment and reagents for LJ DST
- Troubleshooting of non-interpretable LJ DST results.
- Causes of false sensitive and false resistance results in a DST lab.
- Reasons for repeating DST in the laboratory.

# Exercise (5 minutes)

- Joseph a laboratory technologist receives a presumptive MDR sample in the laboratory and tests it on LPA, Gene xpert and 1<sup>st</sup> LJ DST as per the Laboratory's algorithm. The results for all the tests are described in the table below:

GENE XPRT	LPA	1 <sup>ST</sup> line phenotypic DST
MTB DETECTED HIGH, RIFAMPICIN RESISTANCE DETECTED	Isoniazid: RESISTANT Rifampicin: RESISTANT	Pyrazinamide; SENSITIVE Isoniazid; SENSITIVE Rifampicin; SENSITIVE Ethambutol; SENSITIVE

# Exercise (5 minutes)

- Explore all possible causes of such a discordant result and how you would manage such discrepancies in your laboratory?

# The right equipment/materials used correctly:

- Start with the right equipment/supplies
  - Use only sterile tubes and pipettes
  - Sterile glass beads to break up clumps in inoculum from solid media
  - Calibrated pipettes
    - Drug addition to tubes
  - 0.5 McFarland standard with same tube diameter as inoculum tube and same media
    - Organism in 7H9, use 7H9 blank in densitometer
- Vortex well and allow clumps to settle
  - Do not select inoculum from bottom of tube or over inoculate tubes with clumps of organisms



# Use the correct reagents

- Drug reconstitution:
  - Use sterile distilled/deionized water and other diluents.
  - Ensure drug dissolves completely.
  - Store reconstituted drug at  $-20^{\circ}\text{C}$  up to 6 months, but do not exceed original expiration date.
  - Use thawed, reconstituted drug the same day; discard any unused portions.

# Use the correct inoculum

- Start with the right inoculum
  - Pure cultures only
  - Fresh cultures
    - Slant  $\leq$  14 days
  - Homogeneous inoculum
    - Use sterile glass beads and vortex to break up organism clumps
    - Use only supernatant for inoculation,
- Too high or too low an inoculum may give erroneous results or un-interpretable results
  - Prepare dilutions according to procedure.
  - Use accurate pipette to add inoculum to tubes



# DST troubleshooting: Non-interpretable results

## *Contaminated result;*

- Inoculum contained mixed growth of *MTB* contaminants.
- Materials used were not sterile.
- Cross contamination during manipulation.

## *Heavy growth on growth control;*

- Inoculum was heavy.

# LJ DST RESULT; NON INTERPRETABLE RESULTS

- ***No growth on growth control:***
  - Inoculum was too light
  - Inoculated culture isolate was too old containing only dead *Mtb* cells.
- ***Growth on PNB:***
  - Inoculum contained NTM

# LJ DST TROUBLESHOOTING

- **False susceptible**

- Organism inoculum is too light
  - Starting suspension must be 1.0 McFarland
- Suspension not homogeneous
- Culture too old
- Growth control dilution was too heavy
  - Growth control grew too fast
- Was proper amount of bacteria suspension added to the drug containing tubes?
- Drug not reconstituted properly
- A high drug critical concentration in the drug containing media
- Organism could be borderline-resistant

# LJ DST TROUBLESHOOTING

- **False resistant**

- Culture is contaminated or mixed culture
- Organism inoculum was too heavy
- Organism was not properly diluted
- Suspension was not homogeneous or had clumps
- Growth control dilution was too light
  - Growth control grew too slowly
- Drug was not reconstituted properly
- Drug was not stored properly
- Drug was not added to LJ media.
- Organism could be borderline-resistant

# LJ DST Troubleshooting: when to retest

- Un-interpretable results
- Mono-resistance (especially RIF, EMB & PZA, CAP, AMK)
- Results disagree with previous testing from the same patient.
  - Borderline resistance/susceptibility?
- QC organisms fail to give expected results
- Unexpected high number of resistant isolates in a batch.

# Assesment

1. List 5 precautions one should take to ensure that accurate DST results are produced in the lab.
2. A turbid McFarland suspension beyond 1.0 can cause a false resistant LJ DST result True/False
3. List 4 causes of false resistant LJ DST results.
4. List 4 causes of false sensitive LJ DST results.



# Summary

- Organism suspension must be homogeneous, without clumps
- Care must be taken when preparing dilutions
- Drugs must be reconstituted and stored properly
- Sterile tubes and pipettes must be used
- Accurate pipettes must be used, not disposable transfer pipettes
- QC of DST is critical to ensure test is functioning properly
  - Should be performed on all new lots of DST drugs.

# References

- GLI TB training package  
<http://www.stoptb.org/wg/gli/trainingpackages.asp>

# Acknowledgments

