

MGIT CULTURE

Module 8: Overview of liquid culture using the **BACTEC[™] MGIT[™] 960 system**

Date:

Venue:

Presenter:

Outline

- Principle of MGIT Culture
- Indicators for MGIT Culture
- Equipment required
- Reagents in the MGIT tube
- Instrument overview
- Loading the inoculated MGIT tube into the instrument

Principle of the analysis

- The BACTEC™ MGIT™960 liquid culture system is based on fluorescence detection of mycobacterial growth in clinical samples other than blood.
- The MGIT 7ml tubes contain modified middle brook 7H9 broth and a fluorescent sensor at the bottom which responds to the concentration of oxygen.
- (An enrichment, MGIT OADC (Oleic acid, Albumin, Dextrose and Catalase) or MGIT 960 Growth Supplement, is added to make the medium complete. This Growth



....CONT (Principle of the analysis)

- Supplement is essential for growth of many mycobacteria, especially those belonging to *M*
- *Tuberculosis* complex. Addition of the MGIT PANTA is necessary to suppress contamination)
- Initial concentration of dissolved oxygen quenches the emissions from the compound and little fluorescence can be detected. **Later**, (On inoculation with a specimen), actively respiring microorganisms consume the oxygen which allows the compound to fluoresce.

....CONT (Principle of the analysis)

- MGIT tubes are incubated into the BACTEC™ MGIT™960 system which monitors the tubes for increasing fluorescence for positive samples or **6 weeks for negative ones (and flags a tube as negative if no fluorescence is detected within 6 weeks (42 days))**
- Positive samples are usually detected from 4 days onwards.

Indications for MGIT Culture

- Diagnosis of TB
- MDR follow-up

EQUIPMENT REQUIRED

- 1000ul adjustable pipette
- Sterile barrier filter pipette
- Mycobacteria Growth Indicator Tubes.
- BACTEC™ MGIT™960 System
- Biological safety Cabinet
- Light Microscope
- Disposable 50ml falcon tubes
- Centrifuge with a minimum of 3000x g force with safety shield

EQUIPMENT REQUIRED (...cont)

- Vortex mixer, Fridge, water bath
- Timer
- Petri dishes
- Markers
- Incubator (35 - 38°C)
- Microscope slides
- Universal bottles
- Sterile plastic loops
- Slide warmer
- Pencil

Reagents

- BACTEC MGIT960 growth Supplement
- BBL MGIT Tube (7 ml)
- MGIT PANTA Antibiotic Mixture
- ZN reagents
- BA reagents
- Formal Milk

BACTEC MGIT 960 reagents

1. MGIT 960 Growth supplement
2. BBL™ MGIT PANTA
3. BBL MGIT Tube (7 ml)



BD Diagnostic Systems

MGIT 960 Growth supplement

- **Oleic Acid**

- Important in mycobacterial metabolism

- **Albumin** (bovine)

- Binds free fatty acids which may be toxic to TB

- **Dextrose**

- Energy source

- **Catalase**

- Destroys toxic peroxides

- **Polyoxyethylene stearate (POES)**

- Enhances growth of *M. tuberculosis* and assists in providing a uniform inoculum

BBL MGIT PANTA

- Antimicrobial mixture:
 - Polymyxin B - Gram negatives
 - Amphotericin B - Yeast, fungi
 - Naladixic Acid - Gram negatives
 - Trimethoprim - Gram negatives and positives
 - Azlocillin - Gram negatives and positives
- Suppresses non-mycobacterial growth

BBL MGIT tube (7ml)

- Modified Middlebrook 7H9 broth base
- Fluorescent indicator detects O₂ level
 - Embedded in silicon in bottom of tube
- Bacteria use O₂ during growth
 - O₂ released from indicator in tube bottom
 - Fluorescence indicates growth of bacteria




**O₂ bound to
indicator = no
fluorescence**

**O₂ consumed during
bacterial growth
= fluorescence**

Negative Culture

Positive Culture



Inoculated tubes ready to incubate

- Place tubes in rack and carry to instrument for loading
 - To prevent breakage due to dropping tubes, always transport tubes in rack
- Load tubes into instrument on same day as processing/inoculation

BACTEC MGIT 960 system

- Incubates cultures at 37 °C
- Continuous monitoring of fluorescence
- Includes bar-code scanning
- 3 drawers with total of 960 tube capacity



BACTEC MGIT 960 system

- Simple workflow
 - 🌐 Instrument automatically directs placement of each tube into instrument
 - 🌐 Instrument indicates positives/Neg with both a visual and an audible signal as they occur
- Automated quality control of instrument
 - Performed continuously to ensure precise and reliable operation.

BACTEC MGIT 960 system

Step 1: Open the desired drawer and select workflow on LCD display



BACTEC MGIT 960 system

Step 2: Scan pre-affixed bar-code on tube




BACTEC MGIT 960 system

Step 3: Load where indicated by solid green LED light in drawer and close the drawer




BACTEC MGIT 960 system

Step 4: Remove positives and completed negatives
as they occur (icons appear on each
drawer)




Notification of positives on drawer

- Indicator lamp on drawer illuminates
- Audible alert sounds




Notification of positives on LCD display



Detection system

- Detector board moves under tubes once an hour to scan for fluorescence
- Data evaluated and tube positivity established based on software algorithms



Length of incubation

- Protocol length for growth detection can be from 1 to 56 days
 - Default protocol length is 42 days for detection of TB from processed sputum specimens
- Specimens with no growth after 42 days can be reported as negative
 - Subcultures or AFB stains are not required unless flakes of growth are seen in tube

BACTEC MGIT 960

Maintenance

Daily

- Check temperature readout on LCD display
- Check temperature of internal thermometer
- Test all three external indicator lamps on all three drawers
- For each drawer, test red and green light-emitting diode (LED) lights
- Check printer paper supply
- Print out "Quality Control Report."

BACTEC MGIT 960 maintenance

Monthly; more often if dusty environment
Clean and replace filters

As needed

Replace fluorescence calibrators prior to
expiration

Quality Assurance Activities Overview

- Do not use expired product
- Reconstituted PANTA must be used within 5 days
- Do not store MGIT tube after addition of PANTA mixture but should immediately be inoculated with the specimen
- Make sure all additions are made inside the biosafety cabinet.
- Do not open several tubes at a time
- Open MGIT tube for as short a period of time as possible

Recap the tube tightly

... CONT

- Do not add more than 500ul of specimen as this may disturb the pH of the medium and cause false fluorescence
- lot tested reagents before use
- Include routine IQC in processes

Summary

- BACTEC MGIT 960 system provides rapid time to detection for TB growth detection
- 960 is easy to use with little hands on time; frees staff to do other tasks
- 960 has built-in quality control for instrument detectors, resulting in minimal maintenance for staff

Assessment review

- Describe the principle of MGIT Culture?
- When is MGIT Culture requested for?
- List 8 equipment required for MGIT culture?

Acknowledgement

