

Training on Proficiency Testing Scheme GeneXpert DTS

Module 1: Overview of proficiency testing

Venue:

Presenter:

Date:

Introduction

Proficiency testing is a essential part of the quality management system that supports laboratories in checking and maintaining quality.

Objectives

By the end of this module participants should be able to

- Define EQA
- Understand the different proficiency testing schemes implemented in TB network
- Understand the different factors to consider before establishing a PT scheme

Module outline

- Introduction
- Uganda SRL scenarios
- Getting started: Issues to Consider for a PT scheme

The Need for Better Diagnostics

- TB remains one of the world's top infectious killers of our time
- Many people die from TB due to delayed diagnosis and treatment initiation
 - Demonstrates the need for more rapid diagnostics



THE REPUBLIC OF UGANDA

Global tuberculosis report 2018. Geneva: World Health Organization; 2018. WHO/CDS/TB/2018.20.
https://www.who.int/tb/publications/global_report/en/

UGANDA Supranational[®]
Reference Laboratory

Timely Accurate Diagnostics for a TB-Free Africa

Working to End TB

The World Health Organization (WHO) End TB Strategy calls for an end to the global tuberculosis (TB) epidemic, aiming to:

- Reduce deaths by 95%
- Cut new TB cases by 90%,
- Ensure that no family is burdened with catastrophic expenses due to TB



THE REPUBLIC OF UGANDA

End TB Strategy - global strategy and targets for tuberculosis prevention, care and control after 2015.

www.who.int/tb/strategy/End_TB_Strategy.pdf

UGANDA Supranational[®]
Reference Laboratory

Timely Accurate Diagnostics for a TB-Free Africa

The End TB Strategy

WHO-recommended rapid TB diagnostics (WRDs) should be available to all persons with signs or symptoms of TB

- All bacteriologically confirmed TB patients should receive drug-susceptibility testing (DST) at least for rifampicin (RIF)
- All patients with RIF-resistant TB should receive DST at least for fluoroquinolones and second-line injectable drugs



THE REPUBLIC OF UGANDA

Global tuberculosis report 2018. Geneva: World Health Organization; 2018. WHO/CDS/TB/2018.20.

https://www.who.int/tb/publications/global_report/en/

UGANDA Supranational[®]
Reference Laboratory

Timely Accurate Diagnostics for a TB-Free Africa

Cepheid Xpert MTB/RIF

The Xpert MTB/RIF® test is a cartridge-based, automated WRD run on the GeneXpert® platform (Cepheid Inc.; Sunnyvale, CA, USA).

- The test can simultaneously detect *Mycobacterium tuberculosis* complex bacteria (MTBC) and resistance to RIF in less than two hours.



Potential Gains

The Xpert MTB/RIF test has the potential to:

- Significantly decrease diagnostic delay
- Increase the detection of drug resistance
- Impact TB transmission



THE REPUBLIC OF UGANDA

UGANDA Supranational[®]
Reference Laboratory

Timely Accurate Diagnostics for a TB-Free Africa

Significant Challenges

There are a variety of challenges to providing quality Xpert MTB/RIF results:

- Inadequate training and mentoring
- Lack of, or poor adherence to, standard operating procedures (SOPs)
- Stock-outs and use of expired reagents
- Inadequate maintenance of equipment
- Poor quality of samples being tested
- Lack of regular on-site supportive supervision

• Lack of monitoring and evaluation of the TB diagnostic network



Introduction to Quality Assurance and Continuous Quality Improvement

Quality assurance (QA) is a system that monitors the various aspects of a diagnostic process ensuring that the results it produces are *accurate, reliable* and *timely*.

Implementation of quality assurance (QA) activities across the TB diagnostic network is part of the continuous quality improvement (CQI) process. CQI is a cyclical, continuous process-based, data-driven approach to improving the quality of diagnostic testing.

CQI operates under the belief that there is always room for improving operations, processes, and activities to increase quality.



THE REPUBLIC OF UGANDA

UGANDA Supranational[®]
Reference Laboratory

Timely Accurate Diagnostics for a TB-Free Africa

CQI: Two Parts

- **National and Supervisory Levels**
 - Focuses on establishing or integrating Xpert MTB/RIF QA activities into the TB diagnostic network in a country or region.
- **QA at Xpert MTB/RIF Testing Sites**
 - addresses key activities to be carried out at the testing sites to ensure the production of quality Xpert MTB/RIF results

Part 1: National and Supervisory Levels

Pillars of a Quality Assurance System

National Planning

- Governance
- Situational analysis
- Planning & budgeting

Procedures & training

- Quality procedures & documentation
- Training & certification
- Clinical–laboratory interface

Testing infrastructure & supplies

- Safe & functional site
- Equipment & supplies
- Data connectivity

Monitor & evaluate

- Remote monitoring
- External quality assessment (EQA)
- M&E



Xpert MTB/RIF PT Panel Availability

- Availability of affordable proficiency testing (PT) material was identified as an important gap in CQI programs
 - In 2012 less than 43% (387/907) of global Xpert MTB/RIF testing sites participated in PT
 - Only 24% (77/317) of Xpert testing sites in the African region participated in PT

WHO TB Laboratory Diagnostic Services Dataset (<https://www.who.int/tb/country/data/download/en/>)




Overview of Quality Assessment

Assessment is a critical aspect of laboratory quality management, and it can be conducted in several ways with the commonest being EQA;

- **Define EQA??**

Method that allows for comparison of a laboratory's testing to a source outside the laboratory

Types of EQA



What is Proficiency Testing?

- Is an evaluation of participant performance against pre-established criteria by means of inter-laboratory comparison.
- A PT program is a quality assessment tool that provides a retrospective measure of technical quality and verifies the reliability of lab results
- However such programs are expensive to establish hence meagre in developing countries.

Types of PT schemes

- **Quantitative scheme** – Quantifies one or more measurands of the proficiency test item;
- **Qualitative scheme** – Identifies or describes one or more characteristics of the proficiency test item. The results of qualitative tests are descriptive and reported on a categorical or ordinal scale, e.g. identification of micro-organisms, or by identifying the presence of a specific measurand (such as a drug or a grading of a characteristic).
- **Simultaneous scheme** – where proficiency test items are distributed for concurrent testing or measurement within a defined time period
- **Continuous scheme** – where proficiency test items are provided at regular intervals



PT schemes for laboratories

- PT is an EQA program involving sending a panel of sample to a group of participating labs.
- Means of verifying the reliability of lab results
- The main objective is to assess proficiency of the labs in the diagnostic TB sites

Objectives of EQA schemes for laboratories

- Evaluation of Lab testing competence
- Assessment of individual testing performance of the lab staff
- Evaluation of reliability of a test procedure
- Production of relevant information to participating labs when necessary



Objectives of EQA schemes for laboratories

- Providing the basis for the comparability of results during epidemiological surveillance and disease control
- Collecting information on laboratory measurements (intra- and inter-laboratory)
- Collecting information for the purpose of licensing or accreditation of laboratories

Group Exercise

- In a group of 3 or individually list the EQA programmes available in your country and explain the limitation of each. (5 minutes)



Uganda SRL PT scheme overview

- The Uganda National Tuberculosis Reference Laboratory (NTRL) Proficiency testing scheme covers
 - Microscopy,
 - Culture,
 - Line probe Assays,
 - GeneXpert MTB/RIF,
 - and Drug susceptibility testing techniques.
 - SARS-COV-2

Design of the NTRL PT scheme;

- randomly selected sub-samples from a bulk homogeneous supply of material

- distributed simultaneously to participating



Uganda SRL PT schedule

EQA panels provided	Microscopy	LPA	Genexpert	Culture	DST	Sars-COV 2
Frequency/ year	Biannual	Annual	Biannual	Biannual	Annual	Biannual
Sending out month	February and August	August	February and August	February and August	August	February and August
Quantity	10 slides/round	10 Dried Tube Specimens /round	05 Dried Tube Specimens /round	10 isolates/round	20 isolates/round	05 Dried Tube Specimens / round
Expected TAT (from date of Dispatch)	44 days	44days	44 days	91 days	116 days	44 days

Uganda SRL: GeneXpert PT Scheme

- Materials; Standard control strains of *Mycobacterium tuberculosis* (MTB) eg H37Rv , Non-MTB strains (NTM) PBS, Cryovials, Cryo labels, Ziplock bags etc,
- Composition; 5 PT panels
- Must be **non-viable**
- Frequency;
 - Availability of funds and purpose
 - Twice a year

Uganda SRL: LPA PT Scheme

- Materials; Standard control strains of *Mycobacterium tuberculosis* (MTB) eg H37Rv , Non-MTB strains (NTM) PBS, Cryovials, Cryo labels, Ziplock bags etc,
- Composition; 5 PT panels
- Must be **non-viable**
- Frequency;
 - Availability of funds and purpose
 - Twice a year

Uganda SRL: Microscopy PT Scheme

- Materials; Standard control strains of *Mycobacterium tuberculosis* (MTB) eg H37Rv , Non-MTB strains (NTM) PBS, Cryovials, Cryo labels, Ziplock bags etc,
- Composition; 10 PT panels At least 10 slides provides a valid and fair test Batch of stained and unstained smears
- Must be **non-viable**
- Frequency;
 - Availability of funds and purpose
 - Twice a year

Uganda SRL: Culture PT Scheme

- Materials; Standard control strains of *Mycobacterium tuberculosis* (MTB) eg H37Rv , Non-MTB strains (NTM) PBS, Cryovials, Cryo labels, Ziplock bags etc,
- Composition; 10 PT panels
- Must be **viable**
- Frequency;
 - Availability of funds and purpose
 - Twice a year

Uganda SRL: SARS Cov -2 PT Scheme


- Materials; Standard control strains of *SARS-Cov 2* PBS, Cryovials, Cryo labels, Ziplock bags etc,
- Composition; 5 PT panels
- Must be **non-viable**
- Frequency;
Availability of funds and purpose
 - Twice a year

Uganda SRL: DST PT Scheme

- Materials; Standard control pure strains of *MTB*, PBS, Cryovials, Cryo labels, Ziplock bags etc,
- Composition; 20 PT panels
- Must be **viable**
- Frequency;
Availability of funds and purpose
 - Once a year

Getting Started: Issues to Consider (1)

- Capacity of NRL
 - HR
 - Training and expertise
- Type and Composition of PT panels
- Internal QC procedures to ensure
 - Homogeneity
 - Stability
- System for sending PT panels
- Frequency of testing
- Forms to record and report results
- Time allowed for technicians to complete PT
- Performance criteria
- Feedback and corrective action if needed
- Mechanism to resolve discrepant results



Getting Started: Issues to Consider (2)

Dispatch of PT panels

- Delivery system based on services, regulations, resources available:
 - mail/post
 - courier
 - supervisory visit
- Turnaround time; consideration of
 - total no. of days from dispatch of panels to closing dates,
 - preparation of individual and summary reports,
 - presentations to PT team,
 - dispatch of the individual and summary reports

Getting Started: Issues to Consider (3)

- Safe package to prevent breakage/damage of PT materials:
 - strong plastic slide holders
- A standardized PT reporting form / an accompanying letter to provide instructions
- No incentives or punitive actions as a result of the PT exercise
- Time allowed to complete the PT exercise

Getting Started: Issues to Consider (4)

Analysis of PT Results

- A scoring system is to be developed prior to test
- Distinguish major and minor errors
- Determine successful score
- Determine plan of action for poor performances
- Database for analysis

Getting Started: Issues to Consider (5)

Feedback to Laboratories on PT results

- Timely and confidential
- Individual and aggregate test results
- Criteria for acceptable performance
- Reports to TB program coordinator should provide appropriate background information and recommendations and not simply scores
- Poor performance often requires a visit to laboratory

Getting Started: Issues to Consider (6)

Feedback to Laboratories on PT results

- Timely and confidential
- Individual and aggregate test results
- Criteria for acceptable performance
- Reports to TB program coordinator should provide appropriate background information and recommendations and not simply scores
- Poor performance often requires a visit to laboratory

Roles of Participating labs

- . To be successful, PT instructions must be followed carefully, all paper work completed accurately, and results submission deadlines met.
- All PT results, as well as corrective actions, should be recorded and the records maintained for an appropriate period of time.
- There must be no difference in the treatment of PT samples and the patient's sample
- Avoid discussion of results with other laboratories.

Limitations Of PT Schemes

- PT results are affected by variables not related to patient samples Such as preparation of the sample, clerical functions, selection of statistical methods of evaluation.
- PT will not detect all problems in the laboratory, particularly those that address the pre- and post-examination procedures.
- A single unacceptable result does not necessarily indicate that a problem exists in the laboratory.

Assessment

- Define the term EQA and discuss its importance
- What factors might drive the choice of an EQA program to implement
- List the 3 EQA programmes that a TB network can implement

Summary

- EQA plays a crucial role in evaluating the performance of a laboratory network
- Proper planning and considerations are vital to drive the choice of an EQA program to implement



References

- External Quality for AFB Smear microscopy by IUATLD
- ISO 13528:2005, *Statistical methods for use in proficiency testing by interlaboratory comparisons*
- ISO 15189:2012, *Medical laboratories – Particular requirements for quality and competence*
- ISO Guide 34, *General requirements for the competence of reference material producers*
- ISO Guide 35, *Reference materials – General and statistical principles for certification*
- ISO/IEC 17043 First edition 2010-02-01
- Guide 34, ISO Guide 35 and ISO 13528 (homogeneity and stability)
- ISO/IEC Guide 98-3, *Uncertainty of measurement – Part 3: Guide to the expression of uncertainty in measurement* (GUM:1995)
- ISO/IEC 17011:2004, *Conformity assessment – General requirements for accreditation bodies accrediting conformity assessment bodies*
- ISO/IEC 17025, *General requirements for the competence of testing and calibration laboratories*



Acknowledgments

