

Module 1: Over view

Uganda Supranational Reference Laboratory



Course outline



- THE REPUBLIC OF UGANDA

 IN COOLUCTION
- Objectives
- Getting started
 - PT schedule
 - Enrolment
 - Sending
 - Performing
- Results analysis
- ISO 17043 standard





- PT is an EQA programme involving sending a panel of sample to a group of participating labs.
- Means of verifying the reliability of lab results
- The main objective is to assess proficiency of the las in the diagnosing TB
- However such programs are not available in developing countries.





In a group of 3 list the EQA programmes available in your country and explain the limitation of each. (5 minutes)





- Evaluation of Lab testing competence
- Assessment of individual testing performance of the lab staff
- Evaluation of reliability of a test procedure
- Production of relevant information to participating labs when necessary



- Providing the basis for the comparability of results during epidemiological surveillance and disease control
- Collecting information on laboratory measurements (intra- and interlaboratory)
- Collecting information for the purpose of licensing or accreditation of laboratories



Introduction



- The Uganda National Tuberculosis Reference Laboratory (NTRL) Proficiency testing scheme
 - microscopy,
 - culture,
 - line probe Assays,
 - GeneXpert MTB/RIF,
 - and drug susceptibility testing techniques.
- design of the NTRL PT scheme; randomly selected sub-samples from a bulk homogeneous supply of material distributed simultaneously to participating laboratories.



Uganda SRL PT sche



Timely Accounts Diagraphatics for a TD Erro Africa

EQA panels category	Microscopy	LPA	GeneXpert	Phenotypic DST
Frequency/ year	Biannual	Biannual	Biannual	Annual
Sending out Month	February and August	February and August	February and August	August
Quantity	10 slides/round	10 isolates/round	04 isolates/roun d	20 isolates/roun d
Expected TAT	10 days	10 days	10 days	3 months




etting Started: Issues to C



- Capacity of NRL
 - •HR
 - Training and expertise
- Type and Composition of PT panels
- Internal QC procedures to ensure
 - Homogeneity
 - Stability
- System for sending PT panels



- Forms to record and report results
- Time allowed for technicians to complete PT
- Availability of microscopes
- Performance criteria
- Feedback and corrective action if needed
- Mechanism to resolve discrepant results





Sending PT panels



- Delivery system based on services, regulations, resources available:
 - mail/post
 - courier
 - supervisory visit
- Turnaround time
- Safe package to prevent breakage/damage of PT materials:
 - strong plastic slide holders





- Frequency
- A standardized PT reporting form / an accompanying letter to provide instructions

- No incentives or punitive actions as a result of the PT exercise
- Time allowed to complete the PT exercise





- A scoring system is to be developed prior to test
- Distinguish major and minor errors
- Determine successful score
- Determine plan of action for poor performances
- Database for analysis





- Timely and confidential
- Individual and aggregate test results
- Criteria for acceptable performance
- Reports to TB program coordinator should provide appropriate background information and recommendations and not simply scores
- Poor performance often requires a visit to laboratory



GeneXpert network – PT panels and support



- \bullet Materials; Standard control strains of <code>Mycobacterium tuberculosis</code> (MTB) H37Rv , Negative stock Cryovials, FTA cards
- Composition; 4 PT panels for each four module machine
 - Must be non-viable
 - Frequency;
 - Availability of funds and purpose
 - At least once a year



Microscopy PT



- Panels' Composition
- The composition of a panel set is determined by NRL
- Number and types of slides to reassure that correct or incorrect results are not accidental
 - At least 10 slides provides a valid and fair test
 - Batch of stained and unstained smears
- Unstained smears:
 - Evaluate staining technique; provide information about stain preparation and quality




ture and DST labs – PT panels and support



Timely Accurate Diagonostics for a TB-Free Africa

- Started in 2010 with 4 labs within Uganda
- Has expanded to include other countries



Composition of PT materials

C & DST – 20 pure characterized isolates

LPA – 10 isolates

Support for PT provision



Timely Accurate Diagonostics for a TB-Free Africa

Technical requirements

- Personnel, accommodation and environme
- Design of PT schemes
- Choice of method or procedure
- Operation of PT schemes
- Data analysis and evaluation
- Reports
- Communication with participants
- confidentiality



Timely Accurate Diagonostics for a TB-Free Africa

Management requirements

- Organization
- Management system
- Document control
- Review of requests, tenders and contracts
- Sub contracting services
- Purchasing services and supplies
- Complaints and appeals
- Corrective actions
- Preventive actions
- Control of records
- Internal audits
- Management reviews



Assessment



- Define the term EQA and discus its importance
- What factors might drive the choose of EQA programme to implement
- List the PT Scheme that a TB lab can implement



References



- ISO 13528:2005, Statistical methods for use in proficiency testing by interlaboratory comparisons
- ISO 15189, Medical laboratories Particular requirements for quality and competence
- ISO Guide 34, General requirements for the competence of reference material producers
- ISO Guide 35, Reference materials General and statistical principles for certification
- ISO/IEC 17043 First edition 2010-02-01
- Guide 34, ISO Guide 35 and ISO 13528 (homogeneity and stability)
- ISO/IEC Guide 98-3, Uncertainty of measurement Part 3: Guide to the expression of uncertainty in measurement (GUM:1995)
- ISO/IEC 17011:2004, Conformity assessment General requirements for accreditation bodies accrediting conformity assessment bodies
- ISO/IEC 17025, General requirements for the competence of testing and calibration laboratories

Acknowledgments





















